

# STUDY RESTORATION

Shock absorption testing using  
COLTENE BRILLIANT Crios sample crowns

## BRILLIANT Crios

**M. Menini**  
University of Genova, Italy  
Study Report to COLTENE 08/2017

### STUDY AIM

Shock absorbing capacity of BRILLIANT Crios compared to zirconia, silicate ceramic and metal.

### EXPERIMENTAL SETUP

Four BRILLIANT Crios crowns were produced by COLTENE according to the specifications of M. Menini. The experimental setup for testing the crowns was according "Int J Prosthodont 26, 549, Menini et al., Shock Absorption Capacity of Restorative Materials for Dental Implant Protheses: An In Vitro Study". The maximum force transmitted onto a simulated peri-implant bone was determined with following materials:

	BRILLIANT Crios	Ney-Oro CB <sup>3</sup>	Empress 2 <sup>3</sup>	Procera Zirconia <sup>3</sup>
<b>Manufacturer</b>	COLTENE	Dentsply Sirona	Ivoclar Vivadent	Nobel Biocare
<b>Material Class</b>	Composite	Gold Alloy	Silicate Ceramic	Zirconia
<b>E-Modulus /GPa</b>	10.3 <sup>2</sup>	77 <sup>1</sup>	96 <sup>1</sup>	210 <sup>1</sup>

<sup>1</sup> Int J Prosthodont 26, 549, Menini et al., Shock Absorption Capacity of Restorative Materials for Dental Implant Protheses: An In Vitro Study

<sup>2</sup> Menini, report to COLTENE

<sup>3</sup> Not Trademarks of COLTENE

### RESULT

No fractures of the samples occurred during the test. All BRILLIANT Crios crowns showed a shock-absorbing behavior similar to composite resin materials tested in previously published studies. BRILLIANT Crios showed up to 57% lower stress transmission compared to zirconia, around 43% less transmission than silicate ceramic and 19% less compared to metal alloys.

### Maximum Force transmitted to peri-implant bone

