


# UNIVERSAL MIXING TIP

One for all, all for one!



## One new UNIVERSAL MIXING TIP – Always the right choice

The new COLTENE Universal Mixing Tip (shown above) contains a non-glare, see-through casing with white-colored connection plug that replaces four different COLTENE mixing tips (shown below with COLTENE part #).



### Product Benefits

This new system enables dental staff to store and dispense one mixing tip faster and more efficiently.

### USP's

One Universal Mixing Tip – always the right choice!

- Use for all 50ml and 25ml microSystem® cartridges from COLTENE.
- Use for all COLTENE light and regular body, bite registration, and gingival mask materials.
- All common oral and spreader tips of 50ml and 25ml cartridges can be used.
- One Universal Mixing Tip is a practical, flexible, and economical solution.
- Up to 55% less extrusion force is required during application of wash impression material.

### Customer Notice

- DENTAL OFFICE - each affected impression or bite registration pack listed will include an illustrated card explaining the change.



### Order Information

UNIVERSAL MIXING TIP  
40 pcs.

REF 6001 9663

# One new UNIVERSAL MIXING TIP – Always the right choice

## COLTENE Universal Mixing Tip Reference Guide: Refills and Packs

COLTENE Cat. No.	Description	Previous Mix Tip Cat. No.	New Mix Tip Cat. No.	Quantity Mix Tips	Quantity Oral Tips
#4626	President® 50ml LB SGL	6550-YEL	6001 9663	12	0
#4627	President 50ml RB SGL	6550-YEL	6001 9663	12	0
#6002	President microSystem® LB SGL	6202-GRE	6001 9663	12	0
#6012	President microSystem RB SGL	6214-CLEAR	6001 9663	12	0
#6400	Jet Bite 50ml SGL	6551-PINK	6001 9663	12	12*
#6440	Jet Bite microSystem SGL	6214-CLEAR	6001 9663	12	12*
#8060	GI Mask® Automix Starter	6550-YEL	6001 9663	12	12
#8065	GI Mask Automix Refill	6550-YEL	6001 9663	12	12
#6735	Magic Foam Cord® 50ml	6550-YEL	6001 9663	24	24
#6500	Affinis® microSystem LB SGL	6202-GRE	6001 9663	18	18
#6600	Affinis Fast microSystem LB	6202-GRE	6001 9663	18	18
#6610	Affinis Fast microSystem RB SGL	6214-CLEAR	6001 9663	18	18
#6501	Affinis LB Reg Set SGL	6550-YEL	6001 9663	12	0
#6510	Affinis microSystem RB SGL	6214-CLEAR	6001 9663	18	18
#6511	Affinis LB Fast Set SGL	6550-YEL	6001 9663	12	0
#6601	Affinis LB Reg Set SGL	6550-YEL	6001 9663	12	0
#6611	Affinis RB Fast Set SGL	6550-YEL	6001 9663	12	0
#6775	Affinis Precious® LB SGL	6550-YEL	6001 9663	12	0
#6776	Affinis Precious RB SGL	6550-YEL	6001 9663	12	0
#6779	Affinis Precious LB SGL	6202-GRE	6001 9663	18	18
#6780	Affinis Precious RB SGL	6214-CLEAR	6001 9663	18	18
6001 9663	UNIV MIXING TIP REFILL 40 PCS				

\*Qty. Spreader Tips



### Order Information

UNIVERSAL MIXING TIP  
40 pcs.

REF 6001 9663

**Coltene/Whaledent Inc.**  
235 Ascot Parkway  
Cuyahoga Falls, OH 44223 / USA  
Tel. USA & Canada + 1 800 221 3046  
+ 1 330 916 8800  
Fax +1 330 645 8704  
info.us@coltene.com

**Coltene/Whaledent AG**  
Feldwiesenstrasse 20  
9450 Altstätten / Switzerland  
Tel. +41 (0)71 757 53 00  
Fax +41 (0)71 757 53 01  
info.ch@coltene.com



**Coltene/Whaledent GmbH + Co. KG**  
Raiffeisenstrasse 30  
89129 Langenau / Germany  
Tel. +49 (0)7345 805 0  
Fax +49 (0)7345 805 201  
info.de@coltene.com

coltene.com

**COLTENE**